

## Product datasheet for **HK212207**

### CSRP2BP (KAT14) Human qPCR Template Standard (NM\_177926)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_177926
Synonyms:	CRP2 binding partner; CRP2 binding protein; CRP2BP; CSRP2 binding protein; cysteine rich protein 2 binding protein; dj717M23.1; MGC15388; PRO1194
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213718) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)