

## Product datasheet for **HK212252**

### Signal Peptide Peptidase (HM13) Human qPCR Template Standard (NM\_178582)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_178582
Synonyms:	H13; IMP1; IMPAS; IMPAS-1; MSTP086; PSENL3; PSL3; SPP; SPPL1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226073) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)