

## Product datasheet for **HK212264**

### IL22 RA2 (IL22RA2) Human qPCR Template Standard (NM\_181309)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_181309
Synonyms:	CRF2-10; CRF2-S1; CRF2X; IL-22BP; IL-22R-alpha-2; IL-22RA2; ZCYTOR16
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226233) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)