

## Product datasheet for **HK212275**

### CD300E Human qPCR Template Standard (NM\_181449)

#### Product data:

|                  |                                                                                                                                                                                                                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | qPCR Template Standards                                                                                                                                                                                                                                                                       |
| ACCN:            | NM_181449                                                                                                                                                                                                                                                                                     |
| Synonyms:        | CD300LE; CLM-2; CLM2; CMRF35-A5; IREM-2; IREM2; PlgR-2; PlgR2                                                                                                                                                                                                                                 |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (HP218940) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT                                                                                                                                                                                                                      |
| Storage:         | Store at -20°C.                                                                                                                                                                                                                                                                               |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.                                                                                                                                                                                                           |



[View online »](#)