

Product datasheet for **HK212286**

CMTM3 Human qPCR Template Standard (NM_181554)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_181554
Synonyms:	BNAS2; chemokine-like factor super family 3; chemokine-like factor superfamily 3; CKLF-like MARVEL transmembrane domain containing 3; CKLFSF3; FLJ31762; MGC51956; OTTHUMP00000174790
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP229863) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

