

Product datasheet for **HK212394**

TMPRSS4 Human qPCR Template Standard (NM_183247)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_183247
Synonyms:	membrane-type serine protease 2; MT-SP2; MT-SP2, TMPRSS3; TMPRSS3; transmembrane protease, serine 4; transmembrane serine protease 3; type II membrane serine protease
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226543) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)