

Product datasheet for **HK212428**

Nucleoside Diphosphate Kinase 7 (NME7) Human qPCR Template Standard (NM_197972)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_197972
Synonyms:	CFAP67; MN23H7; NDK 7; NDK7; nm23-H7
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210483) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)