

Product datasheet for **HK212459**

PRMT1 Human qPCR Template Standard (NM_198319)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_198319
Synonyms:	6720434D09Rik; ANM1; arginine N-methyltransferase 1; AW214366; HCPI; heterogeneous nuclear ribonucleoproteins methyltransferase-like 2; HRMT1L2; Hrmt1L2; IR1B4; Mrmt1; OTTMUSP00000022387; protein arginine N-methyltransferase 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205398) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

