

## Product datasheet for **HK212479**

### HSD11B1L Human qPCR Template Standard (NM\_198533)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_198533
Synonyms:	11-beta-HSD3; 11-DH3; HSD1L; HSD3; SCDR10; SCDR10B; SDR26C2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP219380) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)