

## Product datasheet for **HK212500**

### **GABA A Receptor gamma 2 (GABRG2) Human qPCR Template Standard (NM\_198903)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_198903
<b>Synonyms:</b>	CAE2; DEE74; ECA2; EIEE74; FEB8; GEFSP3
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP226246) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)