

Product datasheet for **HK212538**

MHF1 (CENPS) Human qPCR Template Standard (NM_199295)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_199295
Synonyms:	apoptosis-inducing, TAF9-like domain 1; CENP-S; CENPS; centromere protein S; MGC32686; OTTHUMP00000065118; OTTHUMP00000065119
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225876) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

