

Product datasheet for **HK212550**

PREI3 (MOB4) Human qPCR Template Standard (NM_199482)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_199482
Synonyms:	2C4D; CGI-95; MOB1; MOB3; MOBKL3; PHOCN; PREI3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP233849) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

