

Product datasheet for **HK212566**

TCF7 Human qPCR Template Standard (NM_201633)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_201633
Synonyms:	FLJ36364; MGC47735; OTTHUMP00000159391; TCF-1; Transcription factor-7, T-cell specific; transcription factor 7 (T-cell specific, HMG-box)
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206769) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)