

Product datasheet for **HK212580**

DPAGT1 Human qPCR Template Standard (NM_203316)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_203316
Synonyms:	GPT, ALG7, DGPT, G1PT, UAGT, UGAT, DPAGT, CDG-lj, DPAGT2, D11S366
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205274) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)