

Product datasheet for **HK212632**

Calpain 3 (CAPN3) Human qPCR Template Standard (NM_212467)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_212467
Synonyms:	CANP3; CANPL3; LGMD2; LGMD2A; MGC4403; MGC10767; MGC11121; MGC14344; nCL-1; p94
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200056) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)