

Product datasheet for **HK212662**

NUDT19 Human qPCR Template Standard (XM_372723)

Product data:

Product Type:	qPCR Template Standards
ACCN:	XM_372723
Synonyms:	androgen regulated gene RP2; D7Rp2; D7Rp2-r; D7Rp2-s; D7Rp2e; nudix (nucleoside diphosphate linked moiety X)-type motif 19; OTTMUSP00000032016; RP2-r; RP2-s
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225860) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)