


Product datasheet for **HK212671**

SARS-CoV-2 (COVID-19) S protein gene qPCR Template Standard

Product data:

Product Type:	qPCR Template Standards
ACCN:	NC_045512.2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP234776) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

Product images:



Relative positions of amplicon targets on SARS-CoV-2 (previous Wuhan-CoV) genome. ORF: open reading frame; RdRp: RNA-dependent RNA polymerase; S: surface glycoprotein; N: nucleocapsid. Numbers below amplicon are genome locations corresponding to Wuhan-CoV genome, MN908947.



[View online »](#)