

Product datasheet for **HK300022**

MIR103-1-as Human qPCR Template Standard (MI0007261)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0007402
Synonyms:	hsa-miR-103-as; MI0007261; MIMAT0007402; MIR103-1-as hsa-mir-103-1-as
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

