

Product datasheet for **HK300033**

MIR1180 Human qPCR Template Standard (MI0006273)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0005825
Synonyms:	hsa-miR-1180; MI0006273; MIMAT0005825; MIR1180 hsa-mir-1180
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)