

Product datasheet for **HK300090**

MIR1266 Human qPCR Template Standard (MI0006403)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0005920
Synonyms:	hsa-miR-1266; MI0006403; MIMAT0005920; MIR1266 hsa-mir-1266
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)