

Product datasheet for **HK300096**

MIR1271 Human qPCR Template Standard (MI0003814)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0005796

Synonyms: hsa-miR-1271; MI0003814; MIMAT0005796; MIR1271 hsa-mir-1271

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)