

Product datasheet for **HK300109**

MIR1282 Human qPCR Template Standard (MI0006429)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0005940

Synonyms: hsa-miR-1282; MI0006429; MIMAT0005940; MIR1282 hsa-mir-1282

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)