

## Product datasheet for **HK300172**

### MIR152 Human qPCR Template Standard (MI0000462)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0000438
Synonyms:	hsa-miR-152; MI0000462; MIMAT0000438; MIR152 hsa-mir-152
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)