

Product datasheet for **HK300176**

MIR1539 Human qPCR Template Standard (MI0007260)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | MIMAT0007401 |
| Synonyms: | hsa-miR-1539; MI0007260; MIMAT0007401; MIR1539 hsa-mir-1539 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)