

## Product datasheet for **HK300190**

### **MIR1827 Human qPCR Template Standard (MI0008195)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | MIMAT0006767  |
| <b>Synonyms:</b>        | hsa-miR-1827; MI0008195; MIMAT0006767; MIR1827 hsa-mir-1827   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| <b>Quality Control:</b> | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment   |
| <b>Storage:</b>         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)