

Product datasheet for **HK300239**

MIR2053 Human qPCR Template Standard (MI0010487)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | MIMAT0009978 |
| Synonyms: | hsa-miR-2053; MI0010487; MIMAT0009978; MIR2053 hsa-mir-2053 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)