

Product datasheet for **HK300293**

MIR302f Human qPCR Template Standard (MI0006418)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0005932
Synonyms:	hsa-miR-302f; MI0006418; MIMAT0005932; MIR302f hsa-mir-302f
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

