

Product datasheet for **HK300296**

MIR30c-1 Human qPCR Template Standard (MI0000736)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0000244
Synonyms:	hsa-miR-30c; MI0000736; MIMAT0000244; MIR30c-1 hsa-mir-30c-1
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)