

Product datasheet for **HK300318**

MIR33B Human qPCR Template Standard (MI0003646)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0003301

Synonyms: hsa-miR-33b; MI0003646; MIMAT0003301; MIR33B hsa-mir-33b

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)