

Product datasheet for **HK300320**

MIR342 Human qPCR Template Standard (MI0000805)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0004694

Synonyms: hsa-miR-342-5p; MI0000805; MIMAT0004694; MIR342 hsa-mir-342

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)