

## Product datasheet for **HK300321**

### MIR345 Human qPCR Template Standard (MI0000825)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | MIMAT0000772  |
| Synonyms:        | hsa-miR-345; MI0000825; MIMAT0000772; MIR345 hsa-mir-345  |
| Components:      | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment   |
| Storage:         | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)