

## Product datasheet for **HK300409**

### MIR515-1 Human qPCR Template Standard (MI0003144)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** MIMAT0002826

**Synonyms:** hsa-miR-515-5p; MI0003144; MIMAT0002826; MIR515-1 hsa-mir-515-1

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)