

## Product datasheet for **HK300413**

### MIR518A1 Human qPCR Template Standard (MI0003170)

#### Product data:

**Product Type:** qPCR Template Standards

**ACCN:** MIMAT0005457

**Synonyms:** hsa-miR-518a-5p; MI0003170; MIMAT0005457; MIR518A1 hsa-mir-518a-1

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)