

Product datasheet for **HK300415**

MIR518c Human qPCR Template Standard (MI0003159)

Product data:

Product Type: qPCR Template Standards

ACCN: MIMAT0002847

Synonyms: hsa-miR-518c; MI0003159; MIMAT0002847; MIR518c hsa-mir-518c

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)