

Product datasheet for **HK300431**

MIR520g Human qPCR Template Standard (MI0003166)

Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0002858
Synonyms:	hsa-miR-520g; MI0003166; MIMAT0002858; MIR520g hsa-mir-520g
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

