

## Product datasheet for **HK300447**

### MIR548A1 Human qPCR Template Standard (MI0003593)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	MIMAT0003251
Synonyms:	hsa-miR-548a-3p; MI0003593; MIMAT0003251; MIR548A1 hsa-mir-548a-1
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)