

Product datasheet for **HK300575**

MIR663B Human qPCR Template Standard (MI0006336)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | MIMAT0005867 |
| Synonyms: | hsa-miR-663b; MI0006336; MIMAT0005867; MIR663B hsa-mir-663b |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)