

Product datasheet for **HK300636**

MIR103-2-as Human qPCR Template Standard (MI0007262)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TCATAGCCCTGTACAATGCTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0007402

Synonyms: hsa-miR-103-as; MI0007262; MIMAT0007402; MIR103-2-as hsa-mir-103-2-as

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)