

## Product datasheet for **HK300640**

### MIR1255B-2 Human qPCR Template Standard (MI0006436)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CGGATGAGCAAAGAAAGTGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0005945
Synonyms:	hsa-miR-1255b; MI0006436; MIMAT0005945; MIR1255B-2 hsa-mir-1255b-2
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)