

Product datasheet for **HK300659**

MIR320c-2 Human qPCR Template Standard (MI0008191)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AAAAGCTGGGTTGAGAGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0005793
Synonyms:	hsa-miR-320c; MI0008191; MIMAT0005793; MIR320c-2 hsa-mir-320c-2
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

