

Product datasheet for **HK300662**

MIR492 Human qPCR Template Standard (MI0003131)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGGACCTGCGGGACAAGAT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0002812

Synonyms: hsa-miR-492; MI0003131; MIMAT0002812; MIR492 hsa-mir-492

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)