

Product datasheet for **HK300701**

MIR1972-2 Human qPCR Template Standard (MI0015977)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCAGGCCAGGCACAGTGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0009447
Synonyms:	hsa-miR-1972; MI0015977; MIMAT0009447; MIR1972-2 hsa-mir-1972-2
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)