

Product datasheet for **HK300729**

MIR3118-6 Human qPCR Template Standard (MI0015981)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGTGACTGCATTATGAAAATTC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0014980

Synonyms: hsa-miR-3118; MI0015981; MIMAT0014980; MIR3118-6 hsa-mir-3118-6

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)