

Product datasheet for **HK300773**

MIR3162 Human qPCR Template Standard (MI0014192)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TTAGGGAGTAGAAGGGTGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0015036
Synonyms:	hsa-miR-3162-5p; MI0014192; MIMAT0015036; MIR3162 hsa-mir-3162
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)