

## Product datasheet for **HK300808**

### MIR3195 Human qPCR Template Standard (MI0014240)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CGCGCCGGGCCCGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0015079
Synonyms:	hsa-miR-3195; MI0014240; MIMAT0015079; MIR3195 hsa-mir-3195
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)