

## Product datasheet for **HK300810**

### MIR3197 Human qPCR Template Standard (MI0014245)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GGAGGCGCAGGCTCGAAAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0015082
Synonyms:	hsa-miR-3197; MI0014245; MIMAT0015082; MIR3197 hsa-mir-3197
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)