

Product datasheet for **HK300834**

MIR3622b Human qPCR Template Standard (MI0016014)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGGCATGGGAGGTCAGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0018005
Synonyms:	hsa-miR-3622b-5p; MI0016014; MIMAT0018005; MIR3622b hsa-mir-3622b
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)