

Product datasheet for **HK300846**

MIR3665 Human qPCR Template Standard (MI0016066)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGCAGGTGCGGGGCGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0018087
Synonyms:	hsa-miR-3665; MI0016066; MIMAT0018087; MIR3665 hsa-mir-3665
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)