

## Product datasheet for **HK300854**

### MIR3674 Human qPCR Template Standard (MI0016075)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	ATTGTAGAACCTAAGATTGGC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0018097
Synonyms:	hsa-miR-3674; MI0016075; MIMAT0018097; MIR3674 hsa-mir-3674
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)