

Product datasheet for **HK300868**

MIR3689d-1 Human qPCR Template Standard (MI0016834)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | GGGAGGTGTGATCTCACACT |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0019008 |
| Synonyms: | hsa-miR-3689d; MI0016834; MIMAT0019008; MIR3689d-1 hsa-mir-3689d-1 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)